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
The screenshot shows the Global Worksheet interface with a plot titled 'Sample-Tube_001'. The plot displays a scatter of points with a red gate labeled 'Parent' at the bottom left. The statistics table below the plot is as follows:

Tube:	Tube_001
Population	#Events %Parent %Total
All Events	10,000 ### 100.0

Annotations in the image include:

- 'Create new global worksheets.' pointing to the Global Worksheets list in the left pane.
- 'Customize plots using the Plot Inspector.' pointing to the Plot Inspector window at the bottom.
- 'Create custom text and graphics.' pointing to the plot area.

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The screenshot shows a plot titled 'Sample-Tube_001' with a red gate labeled 'Parent'. The statistics table below the plot is as follows:

Tube:	Tube_001
Population	#Events %Parent %Total
All Events	10,000 ### 100.0
Parent	1,867 18.7 18.7
Child A	128 8.9 1.3
Child B	219 11.7 2.2

Annotations in the image include:

- 'Use the population hierarchy to verify parent/child relationships.' pointing to the statistics table.
- 'Verify that gates are set appropriately for all samples.' pointing to the gate on the plot.

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The screenshot shows the Batch Analysis dialog box with the following settings:

- Output To Printer:
- Save as PDF:
- Add Report to PDF:
- Statistics:
- Freeze Exponential Scales:
- Use Preferred Global Worksheet:

PDF Filename: sheet\Batch_Analysis_05072007133515.pdf
Export Filename: tistics\Batch_Analysis_05072007133515.csv


Status: 0%

Buttons: Start, Pause, Continue, Close

Annotations in the image include:

- 'Specify where to save the PDF and exported statistics files.' pointing to the filename fields.
- 'Select to print, save as a PDF, or export the statistics as needed.' pointing to the checkboxes.

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BD

Helping all people
live healthy lives

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